# **Specifications**

### Appendix A

#### **Measurement Modes**

% FS Phase sensitive detector output

expressed as a percentage of the present

full-scale sensitivity setting

SIGNAL Phase sensitive detector output

expressed directly in terms of voltage at

input to signal channel

**NOISE** Noise in a bandwidth defined by the output filter

time constant and slope controls and centered at the reference frequency expressed as a percentage of

the present full-scale sensitivity setting

**Harmonic** Fundamental (F) or 2F modes

### **Displays & Indicators**

Two,  $3\frac{1}{2}$  -digit liquid crystal displays, analog center-zero panel meter and back-lit LED indicators show the settings of all the main instrument controls and outputs.

### **Signal Channel**

#### **Voltage Inputs**

Modes A only or Differential (A-B)

Full-scale Sensitivity 100 nV to 3 V rms in a 1-3-10 sequence

Impedance  $100 \text{ M}\Omega \text{ // } 30 \text{ pF}$ 

Maximum Input  $\pm 100 \text{ V DC}$ ; 30 V AC pk-pk without

damage, 10 V AC pk-pk without saturation

Voltage Noise  $5 \text{ nV/}\sqrt{\text{Hz}}$  at 1 kHz typ

CMRR > 100 dB at 1 kHz degrading by 6 dB/octave

Frequency Response 0.5 Hz to 120 kHz

Grounding BNC shields can be grounded or floated via 1 k $\Omega$ 

to ground

#### **Current Input**

Mode  $10^6 \text{ V/A or } 10^8 \text{ V/A}$ 

Full-scale Sensitivity

 $10^{8} \text{ V/A}$  10 fA to 30 nA in a 1-3-10 sequence  $10^{6} \text{ V/A}$  10 fA to 3  $\mu$ A in a 1-3-10 sequence

Frequency Response

10<sup>8</sup> V/A -3 dB at 330 Hz 10<sup>6</sup> V/A -3 dB at 60 kHz

Impedance

 $\begin{array}{ll} 10^8 \text{ V/A} & <2.5 \text{ k}\Omega \text{ at } 100 \text{ Hz} \\ 10^6 \text{ V/A} & <250 \text{ }\Omega \text{ at } 1 \text{ kHz} \end{array}$ 

Maximum Input 15 mA continuous, 1 A momentary

without damage.  $10 \mu A$  AC pk-pk without saturation on  $10^6$  V/A;  $100 \mu A$  AC pk-pk

without saturation on  $10^8 \, \text{V/A}$ 

Noise

10<sup>8</sup> V/A 13 fA/ $\sqrt{\text{Hz}}$  at 500 Hz 10<sup>6</sup> V/A 130 fA/ $\sqrt{\text{Hz}}$  at 1 kHz

Grounding BNC shield can be grounded or floated via 1 k $\Omega$ 

to ground

**Line Notch Filter** > 34dB attenuation  $(a) \pm 1\%$  of 50 or 60 Hz

and/or 100 or 120 Hz

**Dynamic Reserve** 130 dB max

Gain Accuracy

Flat Mode 1% typical Bandpass Mode 2% typical

Gain Stability 200 ppm/°C typical

#### **Reference Channel**

TTL Input

Frequency Range 0.5 Hz to 120 kHz

**Analog Input** 

 $\begin{array}{ll} \text{Impedance} & 1 \text{ M}\Omega \text{ // } 30 \text{ pF} \\ \text{Frequency Range} & 0.5 \text{ Hz to } 120 \text{ kHz} \end{array}$ 

Level

Sinusoidal Input 1.0 V rms\*\*
Squarewave Input 100 mV rms\*\*

\*\*Note: Lower levels can be used with the analog input at the expense of increased phase errors.

Maximum input voltage 5.0 V rms

Phase

Set Resolution 0.1° (front-panel) or 0.005° (computer command

only) increments

Accuracy ±1° typical

Noise 0.005°rms at 100 ms TC, 12 dB/octave

Orthogonality

Above 5Hz  $90^{\circ} \pm 0.5^{\circ}$  0.5Hz - 5Hz  $90^{\circ} \pm 5^{\circ}$  max Drift (Flat Mode)  $< 0.05^{\circ}/^{\circ}$ C

**Lock Acquisition Time** 2 cycles + 100 ms

#### **Demodulator**

**Description** Switching type demodulators operating in either square

wave or Walsh function modes.

**Output Zero Stability** 

High Dynamic Reserve 500 ppm/°C Normal 50 ppm/°C High Stability 5 ppm/°C

**Harmonic Rejection** 

Low-Pass >80 dB at 1 kHz Bandpass >60 dB at 1 kHz

**Time Constant** 

Main output 1 ms to 3 ks in a 1-3-10 sequence

Roll-off 6 and 12 dB/octave
P.S.D. Monitor Output 100 μs nominal
Roll-off 6 dB/octave only

**Offset** Auto and Manual:  $\pm 150 \% FS$ 

Oscillator

Frequency

Range 0.5 Hz to 120 kHz Setting Resolution better than 1%

Absolute Accuracy  $\pm 2\%$ 

**Distortion (THD)** 0.5%

**Amplitude** 

Range

Front panel 1 mV to 1.999 V

Computer Control 1 mV to 2.000 V and 5.000 V

**Setting Resolution** 

1 mV to 500 mV 1 mV 501 mV to 2 V 4 mV

**Output** 

Impedance  $900 \Omega$ 

**Auxiliary Inputs** 

**AUX ADC INPUT CH1 - CH4** 

 $\begin{array}{ll} \text{Maximum Input} & \pm 15 \text{ V} \\ \text{Resolution} & 1 \text{ mV} \end{array}$ 

Input Impedance  $1 \text{ M}\Omega // 30 \text{ pF}$ 

Sample Rate

CH1 only 200 Hz max.

CH1 - CH4 50 Hz max.

Trigger Mode Internal or External Trigger input TTL compatible

#### **Outputs**

**OUT Analog Output** 

Function Output, Noise, Ratio and Log Ratio. Amplitude  $\pm 15 \text{ V } (\pm 10.0 \text{ V} = \pm \text{ full scale})$ 

Impedance  $1 \text{ k}\Omega$ 

**Signal Monitor** 

 $\begin{array}{ll} \text{Amplitude} & \pm 10 \text{ V max} \\ \text{Impedance} & 1 \text{ k}\Omega \end{array}$ 

Aux D/A Outputs

Maximum Output $\pm 15 \text{ V}$ Resolution1 mVOutput Impedance $< 150 \Omega$ 

**Reference Output** 

Waveform 0 to 5 V square wave Impedance TTL compatible

**Power - Low Voltage**  $\pm 15 \text{ V}$  at 100 mA rear panel DIN connector for

powering **SIGNAL RECOVERY** preamplifiers

**Interfaces** 

RS232 and GPIB (IEEE-488). All settings can

be adjusted from the front-panel

General

**Power Requirements** 

Voltage 110/120/220/240 VAC

Frequency 50/60 Hz Power < 130 VA

**Dimensions** 

Width 440 mm (17.25") Depth 89 mm (16.5 ")

Height

With feet 105 mm (4.1 ") Without feet 89 mm (3.5 ")

**Weight** 9.1 kg (20 lbs)

All specifications subject to change without notification

## **Pinouts**

### - Appendix **B**

### **B1 RS232 Connector Pinout**

Figure B-1, RS232 Connector (Female)

Pin	Function	Description
1	Earth Ground	Ties the chassis of the model 5209 to that of the computer
2	Transmit Data	The 5209 transmits data on this line
3	Receive Data	The 5209 receives data on this line
4	Request to Send	This line is asserted by the 5209 when the input buffer is not full
5	Clear to Send	The computer should assert this line to allow the 5209 to transmit data. If left unconnected, the line assumes the asserted state allowing data transmission to proceed
7	Logic Ground	Data signals are referenced with respect to the voltage at this pin

All other pins are not connected

## **B2 Preamplifier Power Connector Pinout**

Figure B-2, Preamplifier Power Connector

Pin	Function
1	-15 V
2	Ground
3	+15 V

Pins 4 and 5 are not connected. Shell is shield ground.